




ATTENTION: Before using SurgyPad, ensure that the OR table mattress is securely affixed to the OR table.



1. Place SurgyPad directly onto the OR table with the black anti-skid strip down & the “Head Side” label up. The SurgyPad must make direct contact with the OR table mattress & should not be used with table covers or sheets.
2. Attach the six small Velcro straps to the surgical bed rail by looping the straps down behind the table rail, back up through the bottom D-ring, & back down towards the rail to affix the ends of the Velcro together.
3. Place the lift sheet on top of the SurgyPad between the Velcro arm protector straps. The lift sheet should be in the small of the patient’s back, allowing the patients skin of the Scapula and Sacrum region to be directly on the SurgyPad.
4. Position the patient on the SurgyPad using the lift sheet to bring them up and off the SurgyPad. Do not drag the patient across the SurgyPad.
5. Tuck the patient’s arms as shown. If the arm protectors are not needed, let them hang off the side or cut them off.
 - a. Place the patient’s arms in between the built-in arm protectors. Wrap the flap closest to the patient’s body first then the outer flap.
 - b. Pull the Velcro straps over and around the arm protectors through the large D-ring.
 - c. Pull the strap back up and affix the ends of the Velcro together. Cut excess Velcro if necessary.
6. If using the body strap, attach as follows:
 - a. Lay the body strap over the patient’s chest with the SurgyPad logo facing upwards.
 - b. Loop the Velcro ends of the body strap through the bed rails and affix the Velcro to the body strap.