


SURGYNEEDLE



The only disposable & clear 5mm aspiration needle!

Aspiration Solution: The SurgyNeedle

The SurgyNeedle is a patent pending aspirating needle used during laparoscopic and robotic procedures. Pulled for each case, but only used when needed, the SurgyNeedle comes sterile packed eliminating dullness, cleaning and sterilization. Our single use product minimizes the risk of contamination in the OR.



Sharp, Beveled Tip

Enables controlled cyst & organ penetration

Semi-Rigid Shaft

Allows for more maneuverability

Clear & Transparent Shaft

Provides greater fluid visualization

Luer-Lock Connection

Provides safe & simple insertion through a trocar

PRODUCT CODES

780-411: 15-gauge aspiration needle

780-412: 21-gauge aspiration needle

Patent Pending • NMSDC Certified • FDA Registered

U.S. Surgitech, Inc. | 470 Mission St. Unit 6, Carol Stream, IL 60188 | Tel: 800.970.4USS | info@ussurgitech.com