


SurgyCut™

NMSDC Certified • FDA Registered

Fully Disposable Suction Electrode



The SurgyCut suction electrode is available in 5mm and either 35cm or 45cm lengths. The 2-in-1 design gives the ability to use the suction feature or the electrode feature without removing any instrument from the trocar.

The SurgyCut suction electrode comes sterile packed and has the ability to screw onto any suction irrigator with the adapter included in each pack. The fully disposable design eliminates the risk of infection and the potential breakdown of insulation.

PRODUCT CODES:

790-103-535: L-Hook with Suction & Electrode 5mm/35cm
790-103-545: L-Hook with Suction & Electrode 5mm/45cm
790-107-535: J-Hook with Suction & Electrode 5mm/35cm
790-107-545: J-Hook with Suction & Electrode 5mm/45cm
790-111-535: L-Wire with Suction & Electrode 5mm/35cm
790-111-545: L-Wire with Suction & Electrode 5mm/45cm
790-115-535: J-Wire with Suction & Electrode 5mm/35cm
790-115-545: J-Wire with Suction & Electrode 5mm/45cm

790-119-535: Curved Blade with Suction & Electrode 5mm/35cm
790-119-545: Curved Blade with Suction & Electrode 5mm/45cm
790-123-535: Straight Blade with Suction & Electrode 5mm/35cm
790-123-545: Straight Blade with Suction & Electrode 5mm/45cm
790-127-535: Spatula with Suction & Electrode 5mm/35cm
790-127-545: Spatula with Suction & Electrode 5mm/45cm



U.S. SURGITECH
Your Source For Medical Products