



U.S. SURGITECH
Your Source For Medical Products

Latex Free • FDA Registered • NMSDC Certified

SurgyMask

The SurgyMask is a safe and easy to apply positioning accessory for facial protection during robotic and laparoscopic procedures in the supine position.

- Lightweight
- Latex free material
- 360 degrees of protection around the face
- Quick on and off mechanism
- Two intubation channels that allow for the endotracheal tube to be placed and secured during the procedure.

The SurgyMask protects from any pressure or weight from surgical instruments or robotic arms that may come in contact with the patient's face during the procedure.



U.S. Surgitech, Inc. | 551 Kimberly Drive, Carol Stream, IL 60188 | Tel: 800-970-4USS | Fax: 800-971-4USS

www.surgymask.com