



U.S. SURGITECH

Your Source For Medical Products

Laparoscopic Instrumentation

800-970-4877

Info@ussurgitech.com

www.ussurgitech.com

U.S. Surgitech, Inc | 551 Kimberly Drive, Carol Stream, IL 60188

Table of Contents

SurgyCut Disposable Laparoscopic Scissor	5
SurgyCut Disposable Laparoscopic Scissor Product Codes	6
SurgyCut Reposable Laparoscopic Scissor	7
SurgyCut Reposable Laparoscopic Scissor Product Codes	8
SurgyCut Laparoscopic Electrode	9
SurgyCut Laparoscopic Electrode Product Codes	10
SurgyCut Laparoscopic Suction Electrode	11
SurgyCut Laparoscopic Suction Electrode Product Codes	12
SurgyNeedle Fully Disposable Aspiration Needle	13
SurgyNeedle Fully Disposable Aspiration Needle Product Codes	14
U.S. Surgitech Contact Information	15



About U.S. Surgitech

U.S. Surgitech began with a revolutionary vision for the surgical consulting industry – to pair excellence in quality and service with unmatched affordability.

Through our unique understanding of product manufacturing and distribution, we strive to provide innovative surgical solutions for each member of our diverse and wide-reaching client base. We distribute a vast selection of products to a varied audience, including outpatient surgery centers, clinics, hospitals, and government institutions – many of whom are at the vanguard of healthcare services with their own reputations to uphold. Our products help make upholding excellence simple.

In an industry that is constantly shifting, adaptability to new best practices and regulations is vital to success. Our cutting-edge services, coupled with our commitment to providing honest solutions, help clients navigate industry changes expertly and quickly. We devote ourselves to the quality of materials and fine craftsmanship of our product, so clients can conduct healthcare services with confidence and trust.

We are experts at providing the excellent quality you need at prices you can afford.



U.S. SURGITECH
Your Source For Medical Products

SurgyCut

Fully Disposable Laparoscopic Scissor

35cm & 45cm Shaft—Maximizing instrument reach depending on the procedure

Delicate Jaw—Allowing for precise dissecting in small pockets

Rotation Knob— Easy and smooth rotation with the index finger

Ergonomic Handle Design— Providing a comfortable grip when handling and cutting

Regular Metz

Mini Metz

Hook

Our fully disposable scissor is **the complete, low maintenance solution**, coming sterile-packed & ready to use, eliminating the need for reprocessing, reducing the risk of dull blades and infection. Engineered to provide a solid, substantial feel in your hands.



SurgyCut

Fully Disposable Laparoscopic Scissor

NMSDC Certified | FDA Registered


Product Codes	
Fully Disposable Standard Metz 35cm	790-201-535
Fully Disposable Standard Metz 45cm	790-201-545
Fully Disposable Mini Metz 35cm	790-203-535
Fully Disposable Mini Metz 45cm	790-203-545
Fully Disposable Hook 35cm	790-205-535
Fully Disposable Hook 45cm	790-205-545

www.ussurgitech.com/surgycut-scissors



SurgyCut

Reusable Laparoscopic Scissor



Rotation Knob

Easy and smooth rotation with the index finger

Ergonomic Handle Design

Providing a comfortable grip when handling and cutting



Regular Metz



Mini Metz



Hook

The SurgyCut reusable scissor provides the perfect solution for reducing waste, being configured as a reusable handle and disposable shaft. Each of the three available SurgyCut scissor shafts connect seamlessly to our durable handle using a ball and socket connection, delivering a peak performance each time.



SurgyCut

Reusable Laparoscopic Scissor

NMSDC Certified | FDA Registered

Product Codes	
Reusable Standard Metz 35cm	790-200-535
Reusable Standard Metz 45cm	790-200-545
Reusable Mini Metz 35cm	790-202-535
Reusable Mini Metz 45cm	790-202-545
Reusable Hook 35cm	790-204-535
Reusable Hook 45cm	790-204-545
Reusable Handle	790-206

www.ussurgitech.com/surgycut-scissors



SurgyCut

Disposable Laparoscopic Electrodes

The SurgyCut disposable electrodes are available in 5mm and either 35 cm or 45 cm lengths. The handle design provides a comfortable grip when used in any position and the PTFE-coated, non-stick tips prevent eschar build-up. All electrodes are fully disposable, eliminating the risk of infection and the potential breakdown of insulation.

Available in
3 Control
Styles

Foot Control

Hand Control

Handle



J-Wire



L-Wire



J-Hook



L-Hook



Spatula

Available in
5 Tips



SurgyCut

Disposable Laparoscopic Electrodes

NMSDC Certified | FDA Registered


Product	Code	Product	Code
L-Hook HC 5mm/35cm	790-100-535	Curved Blade HC 5mm/35cm	790-116-535
L-Hook HC 5mm/45cm	790-100-545	Curved Blade HC 5mm/45cm	790-116-545
L-Hook FC 5mm/35cm	790-101-535	Curved Blade FC 5mm/35cm	790-117-535
L-Hook FC 5mm/45cm	790-101-545	Curved Blade FC 5mm/45cm	790-117-545
L-Hook Handle 5mm/35cm	790-102-535	Curved Blade Handle 5mm/35cm	790-118-535
L-Hook Handle 5mm/45cm	790-102-545	Curved Blade Handle 5mm/45cm	790-118-545
J-Hook HC 5mm/35cm	790-104-535	Straight Blade HC 5mm/35cm	790-120-535
J-Hook HC 5mm/45cm	790-104-545	Straight Blade HC 5mm/45cm	790-120-545
J-Hook FC 5mm/35cm	790-105-535	Straight Blade FC 5mm/35cm	790-121-535
J-Hook FC 5mm/45cm	790-105-545	Straight Blade FC 5mm/45cm	790-121-545
J-Hook Handle 5mm/35cm	790-106-535	Straight Blade Handle 5mm/35cm	790-122-535
J-Hook Handle 5mm/45cm	790-106-545	Straight Blade Handle 5mm/45cm	790-122-545
L-Wire HC 5mm/35cm	790-108-535	Spatula HC 5mm/35cm	790-124-535
L-Wire HC 5mm/45cm	790-108-545	Spatula HC 5mm/45cm	790-124-545
L-Wire FC 5mm/35cm	790-109-535	Spatula FC 5mm/35cm	790-125-535
L-Wire FC 5mm/45cm	790-109-545	Spatula FC 5mm/45cm	790-125-545
L-Wire Handle 5mm/35cm	790-110-535	Spatula Handle 5mm/35cm	790-126-535
L-Wire Handle 5mm/45cm	790-110-545	Spatula Handle 5mm/45cm	790-126-545
J-Wire HC 5mm/35cm	790-112-535	<p>HC = Hand Control</p> <p>FC = Foot Control</p>	
J-Wire HC 5mm/45cm	790-112-545		
J-Wire FC 5mm/35cm	790-113-535		
J-Wire FC 5mm/45cm	790-113-545		
J-Wire Handle 5mm/35cm	790-114-535		
J-Wire Handle 5mm/45cm	790-114-545		




www.ussurgitech.com/surgycut-electrodes

SurgyCut

Fully Disposable Suction Electrode



Open to Use Suction



Close to Expose Electrode

The SurgyCut suction electrode comes sterile packed and has the ability to screw onto any suction irrigator with the adapter included in each pack. The fully disposable design eliminates the risk of infection and the potential breakdown of insulation.

The SurgyCut suction electrode is available in 5mm and either 35cm or 45cm lengths. The 2-in-1 design gives the ability to use the suction feature or the electrode feature without removing any instrument from the trocar.



SurgyCut

Fully Disposable Suction Electrode

NMSDC Certified | FDA Registered

Product Codes	
L-Hook with Suction & Electrode 5mm/35cm	790-103-535
L-Hook with Suction & Electrode 5mm/45cm	790-103-545
J-Hook with Suction & Electrode 5mm/35cm	790-107-535
J-Hook with Suction & Electrode 5mm/45cm	790-107-545
L-Wire with Suction & Electrode 5mm/35cm	790-111-535
L-Wire with Suction & Electrode 5mm/45cm	790-111-545
J-Wire with Suction & Electrode 5mm/35cm	790-115-535
J-Wire with Suction & Electrode 5mm/45cm	790-115-545
Curved Blade with Suction & Electrode 5mm/35cm	790-119-535
Curved Blade with Suction & Electrode 5mm/45cm	790-119-545
Straight Blade with Suction & Electrode 5mm/35cm	790-123-535
Straight Blade with Suction & Electrode 5mm/45cm	790-123-545
Spatula with Suction & Electrode 5mm/35cm	790-127-535
Spatula with Suction & Electrode 5mm/45cm	790-127-545


www.ussurgitech.com/surgycut-electrodes



SurgyNeedle

Fully Disposable Aspiration Needle

The SurgyNeedle is a patent-pending aspirating needle used during laparoscopic and robotic procedures. Pulled for each case, but only used when needed, the SurgyNeedle comes sterile packed always providing a sharp needle, while also eliminating cleaning and sterilization needs. Our single-use product minimizes the risk of contamination and infection during surgery.



Luer-Lock Connection

Provides safe & simple insertion through a trocar

Semi-Rigid Shaft

Allows for more maneuverability

Clear & Transparent Shaft

Provides greater fluid visualization

Sharp, Beveled Tip

Enables controlled cyst & organ penetration

**The only disposable
& clear 5mm
aspiration needle
on the market!**



SurgyNeedle

Fully Disposable Aspiration Needle

Patent Pending | FDA Registered | NMSDC Certified

Product Codes	
SurgyNeedle 5mm/ 15 Gauge/ 40cm	780-411
SurgyNeedle 5mm/ 18 Gauge/ 40cm	780-413
SurgyNeedle 5mm/ 21 Gauge/ 40cm	780-412
SurgyNeedle 5mm/ 23 Gauge/ 40cm	780-414
SurgyNeedle 5mm/ 23 Gauge/ 15cm	780-415

www.ussurgitech.com/surgyneedle



Contact Us

For General Inquiries or to Request a Free Trial

info@ussurgitech.com

To Submit a Purchase Order

sales@ussurgitech.com

For CE Course Inquiries

education@ussurgitech.com

To Contact Us Directly

Phone: 800-970-4877

Fax: 800-971-4877

www.ussurgitech.com

U.S. Surgitech, Inc | 551 Kimberly Drive, Carol Stream, IL 60188





U.S. SURGITECH

Your Source For Medical Products



U.S. Surgitech, Inc | 551 Kimberly Drive, Carol Stream, IL 60188