

Universal SurgeryPad

Item # 780SP1

A better alternative for Robotic & Laparoscopic Surgery


PROCEDURES

Laparoscopic | Robotic | Gynecology | Urology | Colorectal | General | Thoracic | Cardiac

POSITIONING


Trendelenburg | Reverse Trendelenburg | Lateral Tilt

☒ SAFE & SECURE




- Single-use eliminates cross-contamination risk
- High density, latex-free & moisture wicking foam
- Unique integrated arm protectors wrap around patient's arms for protection of arms, hands, and fingers
- Pressure-relieving foam eliminates risk of brachial plexus injury
- Anti-skid strip on the back of the SurgeryPad adds traction between SurgeryPad and mattress
- Innovative 6-strap system further secures SurgeryPad to table

☒ GREATER FLEXIBILITY



- Adjustable arm protectors accommodate a variety of procedures and patient body types.
- Wide Velcro straps secure the arms close to the patient's body
- Adjustable body strap with Y-configuration to accommodate rail accessories/barriers
- 45" length to fit most patient body types

☒ Cost & Time Effective



- Quick and easy setup
- Reduces "never-events"
- Saves facilities 20% to 40% in costs

